

NOS Assay Kit

Catalog No. 781001

TABLE OF CONTENTS

Contents of the Kit	1
Precautions	2
Warranty and Limitation of Remedy	2
If You Have Problems	2
Storage and Stability	2
Additional Items Required	2
About this Assay	3
Protocol: Measurement of Nitric Oxide Synthase Activity	4
Appendix: Stability of Radiolabeled Arginine	6
Preparation of Media and Reagents	7
References	7
Related Products	7
Quick-Reference Protocol	8
Notes	8

CONTENTS OF THE KIT

Number	Item	Quantity
1	Cerebellum extract, rat brain ^a	1 vial
2	Calmodulin	1 vial
3	Reaction buffer ^b	1 vial
4	N ^G -nitro-L-arginine methyl ester HCl	1 vial
5	Homogenization buffer ^b	1 vial
6	Stop buffer ^b	1 vial
7	Equilibrated resin	1 vial
8	Elution buffer ^b	1 plate
9	Spin cups and cup holders ^c	1 cover

^aThis positive control is a homogenate of rat brain tissue resuspended in homogenization buffer.

^bSee Preparation of Media and Reagents (page 6)

^cSufficient spin cups and cup holders are provided for 50 total reactions.

If any of the items listed above are damaged or missing, please contact our Customer Service Department at (800) 364-9897. We cannot accept returns without prior authorization.

PRECAUTIONS

WARNING: THIS PRODUCT IS NOT INTENDED OR APPROVED FOR USE IN HUMANS OR ANIMALS. USE OF THIS PRODUCT FOR HUMAN OR ANIMAL TESTING IS EXTREMELY HAZARDOUS AND MAY RESULT IN DISEASE, SEVERE INJURY, OR DEATH.

- Please read these instructions carefully before beginning this assay.
- For research use only. Not for human or diagnostic use.

WARRANTY AND LIMITATION OF REMEDY

Cayman Chemical Company makes no warranty of any kind, expressed or implied, including, but not limited to, the warranties of fitness for a particular purpose and merchantability, which extends beyond the description of the chemicals on the face hereof, except that the material will meet our specifications at the time of delivery. Buyer's exclusive remedy and Cayman Chemical Company's sole liability hereunder shall be limited to refund of the purchase price of, or at Cayman Chemical Company's option, the replacement of, all material that does not meet our specifications. Cayman Chemical Company shall not be liable otherwise or for incidental or consequential damages, including, but not limited to, the costs of handling. Said refund or replacement is conditioned on Buyer giving written notice to Cayman Chemical Company within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within said thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

IF YOU HAVE PROBLEMS

Our technical support staff may be reached by phone (800-364-9897, 734-971-3335), fax (734-971-3640) or E-Mail (techserv@caymanchem.com) Monday through Friday 8:00 AM to 6:00 PM EST. In order for our staff to assist you quickly and efficiently, please be ready to supply the lot number of the kit (found on the outside of the box).

STORAGE AND STABILITY

This kit will perform as specified if stored as listed in the Contents of the Kit (see above) and used before the expiration date indicated on the outside of the box.

ADDITIONAL ITEMS REQUIRED

1. [^3H]Arginine monohydrochloride [40-70 Ci/mmol, 1 $\mu\text{Ci}/\mu\text{l}$ (Amersham, Arlington Heights, Illinois, Catalog No. TRK698)] -or- [^{14}C]arginine monohydrochloride [>300 mCi/mmol, 50 $\mu\text{Ci}/\text{ml}$ (Amersham, Arlington Heights, Illinois, Catalog No. CFB63)].
2. Reduced nicotinamide adenine dinucleotide phosphate (NADPH*) (Sigma, Catalog No. 1630 or N7505).
3. 10 mM Tris-HCl (pH 7.4).
4. Scintillation fluid and vials.

*NADPH is not stable in solution; therefore, it is not included in the reaction buffer.

ABOUT THIS ASSAY

Our NOS assay kit is a simple, sensitive and specific assay for nitric oxide synthase (NOS) activity. The NOS assay kit is based on the biochemical conversion of L-arginine to L-citrulline by NOS.¹⁻⁶ This reaction, which represents a novel enzymatic process, involves a five-electron oxidation of a guanidino nitrogen of L-arginine to nitric oxide (NO), together with the stoichiometric production of L-citrulline (see Figure 1, page 3). The reaction consumes 1.5 equivalents of reduced nicotinamide adenine dinucleotide phosphate (NADPH)⁷ and also requires molecular oxygen, calcium, calmodulin, and tetrahydrobiopterin (BH_4).⁸

Measuring NOS activity by monitoring the conversion of arginine to citrulline is currently the standard assay for NOS activity in both crude and purified enzyme preparation. Advantages of the NOS assay kit include the use of radioactive substrates ([^3H]arginine or [^{14}C]arginine) that enable sensitivity to the picomole level, as well as the specificity of the assay for the NOS pathway due to the direct enzymatic conversion of arginine to citrulline in eukaryotic cells. Additionally, the easy separation of neutrally charged citrulline from positively charged arginine allows multiple assays to be performed easily.

For routine assays, radioactive arginine is added to protein extracts or purified NOS samples. After incubation, the reactions are stopped with a buffer containing ethylenediaminetetraacetic acid (EDTA), which chelates the calcium required by nNOS and eNOS and, consequently, inactivates the enzyme. In the case of iNOS, the low pH of the stop buffer (pH 5.5) stops the enzyme-catalyzed reaction. Equilibrated resin, which binds to the arginine, is added to the sample reactions and the reactions are then pipetted into spin cups. The citrulline, being ionically neutral at pH 5.5, flows through the cups completely. The NOS activity is then quantitated by counting the radioactivity in the eluate.




Figure 1. NOS catalyzes a 5-electron oxidation of a guanidino nitrogen of L-arginine to generate NO and L-citrulline. L-Hydroxyarginine is formed as an intermediate that is tightly bound to the enzyme. Both steps in the reaction are dependent on calcium and calmodulin.

PROTOCOL: MEASUREMENT OF NITRIC OXIDE SYNTHASE ACTIVITY

Preparation of Extracts from Tissues and Cultured Cells

The citrulline assay has been used to quantitate levels of NOS activity in tissue homogenates from numerous sources including blood vessels, immune cells, visceral organs, nervous tissue, and cultured cells. NOS activity is relatively unstable; therefore, tissues should be harvested quickly after animal euthanasia. If enzyme assays are to be conducted at a later time, it is best to freeze intact tissues or harvested cultured cells prior to homogenization. Wrap the tissues in aluminum foil, flash freeze the tissues in liquid nitrogen, and then store at -80°C .

Extraction of Proteins from Tissues

Note: Because the level of NOS activity will vary between tissues, the volume of 1X homogenization buffer used may require optimization.

1. Prepare an appropriate volume of 1X homogenization buffer (i.e., a 1:10 dilution of the 10X homogenization buffer). Add 10-20 volumes [volume of buffer (ml): weight of tissue (g)] of ice-cold 1X homogenization buffer to a tissue sample (for tissues with low NOS activity, start with 10 volumes of 1X homogenization buffer).
2. Homogenize the tissue using a tissue grinder or an equivalent tissue homogenizer. Keep the tissue homogenate on ice.
3. Pipet 1 ml aliquots of the tissue homogenate into microcentrifuge tubes and spin the tubes in a microcentrifuge at full speed for 5 minutes at 4°C .
4. Transfer the supernatant to fresh microcentrifuge tubes and keep the tubes on ice until use. This extract will be used in the NOS activity assay (see step 3 of Preparing the Reactions, page 4)

Subcellular Tissue Distribution

The subcellular distribution of NOS is tightly regulated in tissues. Endothelial NOS (eNOS or NOS-III) is largely membrane associated as a result of N-terminal myristoylation.^{9,10} Neuronal NOS (nNOS or NOS-I) is largely soluble in adult rat brain, yet in skeletal muscle, it is predominantly associated with membrane fractions.¹¹ The mechanism for membrane attachment of nNOS remains unclear. Inducible NOS (iNOS) is a soluble enzyme. NOS activity in soluble and membrane-associated fractions can be separated by centrifuging the homogenized tissues at $100,000 \times g$ for 60 minutes. The supernatant contains soluble NOS, while the pellet, which is resuspended in homogenization buffer, contains membrane-associated NOS.

Extraction of Proteins from Tissue Culture Cells

Certain cultured cells, such as endothelial cells and activated macrophages, contain NOS, which can be measured using the citrulline assay. The proteins must first be extracted from the cells as follows:

1. Remove the culture media from the tissue culture cells.
2. Wash the tissue culture cells once with phosphate-buffered saline (PBS) and harvest the tissue culture cells in PBS containing 1 mM EDTA. Transfer the tissue culture cells to microcentrifuge tubes.
3. Spin the microcentrifuge tubes in a microcentrifuge at full speed for 2 minutes to pellet the cells.
4. Remove the supernatant from the pelleted cells by vacuum aspiration and then add 100-500 μl of the 1X homogenization buffer to each microcentrifuge tube of pelleted cells. Sonicate briefly to disrupt the cells.
5. Spin the microcentrifuge tubes in a microcentrifuge at full speed for 5 minutes.
6. Separate the supernatant from the pellet and adjust the resulting protein sample to a concentration of 5-10 mg/ml.

Measurement of Nitric Oxide Synthase Activity in Enzyme Extracts

Note: Refer to Appendix: Stability of Radiolabeled Arginine (page 5) for notes on the purification and stability of radiolabeled arginine.

Incubation of the citrulline assay reaction may be carried out for 10-60 minutes at 22-37°C depending on the tissue being used. High levels of nNOS in nervous tissues¹² and skeletal muscle¹¹ permit brief assays (10-15 minutes) of NOS with room temperature incubations. Lower levels of eNOS in vascular tissues require that assays be performed for prolonged periods (60 minutes).

eNOS and nNOS require calcium for enzyme activity; therefore, it is essential to add calcium to experimental assays. A final free calcium concentration of 75 μM is required for optimal NOS activity.¹² When testing NOS activity from tissue extracts, addition of calmodulin to the reaction is not required. However, when testing purified NOS, the addition of calmodulin is required for nNOS and eNOS (for review, see Reference 1, page 6). Crude extracts of macrophage iNOS, such as the 100,000 x g supernatant solution, require Mg^{2+} for maximum activity, while purified iNOS shows no such dependence.

NOS activity in the citrulline assay is defined as counts per minute (cpm) in an incubated test sample as compared to an appropriate blank. The following control reactions can serve as a blank: a reaction that includes 1 mM N^G -nitro-L-arginine methyl ester HCl (L-NAME, a competitive NOS inhibitor provided at a concentration of 10 mM), a reaction in which the extract is boiled prior to the assay, a reaction in which either NADPH or calcium (for nNOS and eNOS) is omitted, or a reaction that is incubated on ice. As for any quantitative enzyme assay, it is important to verify reaction conditions are such that the assay is linear with respect to time and tissue concentration. The NOS in the rat cerebellum extract provided in this kit is linear for at least a 30-minute reaction. Specific activity and substrate affinity of NOS can be assessed by carrying out replicate reactions in the presence of varying amounts of unlabeled arginine. The K_m (Michaelis constant) of NOS is in the range of 2-20 μM .^{7,10,13} Appropriate concentrations of arginine for kinetic studies are 0.1-100 μM .

Preparing the Reactions

Note: Prepare a control reaction with the provided rat cerebellum extract at step 3 below. Thaw the rat cerebellum extract just before using and keep the thawed extract on ice. Do not refreeze the thawed extract. For crude iNOS samples, include magnesium diacetate at a final concentration of 1 mM.

1. Prepare the reaction mixture on ice by adding the following components to a microcentrifuge tube:

Note: The volumes given here yield sufficient reaction mixture for 10 reactions. The reaction mixture can be stored on ice for up to 24 hours.

Inducible NOS is Ca^{2+} independent.

250 μl of 2x reaction buffer (see Preparation of Media and Reagents, page 6)

50 μl of 10 mM NADPH [freshly prepared in 10 mM Tris-HCl (pH 7.4)]

10 μl of [^3H]arginine (1 $\mu\text{Ci}/\mu\text{l}$) or [^{14}C]arginine (50 $\mu\text{Ci}/\mu\text{l}$)

50 μl of 6 mM CaCl_2

40 μl of dH_2O

2. Store the reaction mixture on ice.

Note: To chemically inhibit the control reaction, add 5 μl of L-NAME to the reaction mixture before adding the tissue extract or NOS enzyme preparation.

3. Combine 40 μl of the reaction mixture with 1-10 μl of tissue extract for the experimental reaction. Combine 40 μl of the reaction mixture with 5 μl of the provided rat cerebellum extract for the control reaction.

Note: Recombinant NOS may also be used for a control reaction. (Calmodulin must be added to reaction using recombinant NOS.)

nNOS and eNOS require 0.1 μM calmodulin when assayed as purified enzymes. The addition of calmodulin to tissue extracts is not necessary, but it is recommended. If the addition of calmodulin is required, add calmodulin to a final concentration of 0.1 μM (i.e., add 5 μl to a 50 μl reaction).

4. Incubate the reaction samples at 22-37°C for 10-60 minutes. (For initial experiments, the reaction should be allowed to proceed at room temperature for 30 minutes.)
5. Stop the reaction by adding 400 μl of stop buffer (see Preparation of Media and Reagents, page 6) to the reaction sample.

Processing the Samples

1. Thoroughly resuspend the equilibrated resin provided. Pipet 100 μ l of the equilibrated resin into each reaction sample.
2. Transfer the reaction samples to spin cups and place the spin cups into cup holders.
3. Centrifuge the spin cups and holders in a microcentrifuge at full speed for 30 seconds.
4. Remove the spin cups from the cup holders and transfer the eluate (i.e., the "flowthrough") to scintillation vials. Add scintillation fluid to the vials and quantitate the radioactivity in a liquid scintillation counter.
5. If determining the ratio of unreacted arginine to citrulline is desired, place the spin cups in fresh microcentrifuge tubes and add 400 μ l of elution buffer (see Preparation of Media and Reagents, page 6) to the spin cup.
6. Spin the microcentrifuge tubes (with the spin cups in them) in a microcentrifuge at full speed for 30 seconds.
7. Remove the spin cups and transfer the eluate to scintillation vials. Add scintillation fluid to the vials and quantitate the radioactivity in a liquid scintillation counter.

In the tissues assayed according to these protocols, citrulline is the major radiolabeled compound in the eluate. This can be readily verified using thin-layer chromatography (TLC). Separation of arginine and relevant metabolites is achieved using silica-gel chromatography plates developed with $\text{CH}_3\text{OH}:\text{NH}_4\text{OH}$ (6:1). With this solvent system, arginine migrates at an R_f (the ratio of the distance traveled by a compound to the distance traveled by the solvent) of ~ 0.1 while citrulline migrates at an R_f of ~ 0.5 .

APPENDIX: STABILITY OF RADIOLABELED ARGININE

Radiolabeled Arginine

Prior to initiating the enzyme assays, it is essential to verify the purity of the radiolabeled arginine, otherwise a high blank value for the liquid scintillation counting will greatly reduce the sensitivity of the assay. To assess the blank value, a reaction mixture is applied to the equilibrated resin. Nonadherent radioactivity is eluted with stop buffer, the eluate is collected and the radioactivity is quantitated in a liquid scintillation counter.

1. Prepare a reaction mix by combining the following components in a microcentrifuge tube.
 - 20 μ l of 2X reaction buffer
 - 4 μ l of 10 mM NADPH [freshly prepared in 10 mM Tris (pH 7.4)]
 - 1-10 μ l of [^3H]arginine (1 $\mu\text{Ci}/\mu\text{l}$) or [^{14}C]arginine (50 $\mu\text{Ci}/\mu\text{l}$)
 - dH_2O to bring the total volume to 40 μ l
2. Store the reaction mix on ice
3. Prepare a reaction sample by combining 100 μ l of well-resuspended equilibrated resin and 10 μ l of the reaction cocktail in a microcentrifuge tube.
4. Add 400 μ l of stop buffer to the reaction sample, mix the contents of the reaction sample, and transfer the reaction sample to a spin cup. Place the spin cup in a spin cup holder.
5. Centrifuge the spin cup and holder in a microcentrifuge at full speed for 30 seconds.
6. Collect 100 μ l of the eluate from the spin cup and quantitate the radioactivity in a liquid scintillation counter.

Greater than 95% of the applied radioactivity should be retained by the spin cup. This represents a relatively low blank value. If more than 5% of the radioactivity flows through the spin cup, it is important to purify the arginine prior to conducting the assay. [^3H]Arginine is prone to radiolytic decay and must be purified every 2 months, while [^{14}C]arginine is more stable but much more expensive.

Purification of Radiolabeled Arginine

Radioactive arginine can be purified with the equilibrated resin included in the NOS assay kit as follows:¹⁴

1. Apply the radioactive arginine to 0.5 ml of equilibrated resin in a disposable spin column [e.g., a Poly-Prep[®] chromatography column (Bio-Rad, Richmond, California, Catalog No. 731-1550)].
2. Wash the column with 5 ml of distilled water.
3. Elute the arginine with two 2-ml washes of elution buffer.
4. Lyophilize the arginine and resuspend the arginine in 2% (v/v) ethanol.

PREPARATION OF MEDIA AND REAGENTS

Elution Buffer

1 M NH_4Cl

2x Reaction Buffer

50 mM Tris-HCl (pH 7.4)

6 μM tetrahydrobiopterin (BH_4)

2 μM flavin adenine dinucleotide

2 μM flavin adenine mononucleotide

10x Homogenization Buffer

250 mM Tris-HCl (pH 7.4)

10 mM EDTA

10 mM ethyleneglycol-*bis*(β -aminoethylether)-N,N,N',N'-tetraacetic acid (EGTA)

Stop Buffer

50 mM N-2-hydroxyethylpiperazine-N'-2-ethanesulfonic acid (HEPES) (pH 5.5)

5 mM EDTA

REFERENCES

1. Bredt, D.S. and Snyder, S.H. Nitric oxide: A physiologic messenger molecule. *Annu. Rev. Biochem.* 63, 175-195 (1994).
2. Ignarro L.J. Signal transduction mechanisms involving nitric oxide. *Biochem. Pharmacol.* 41, 485-90 (1991).
3. Marletta, M.A. Nitric oxide synthase structure and mechanism. *J. Biol. Chem.* 268, 12231-12234 (1993).
4. Moncada, S. and Higgs, A. The L-arginine-nitric oxide pathway. *N. Engl. J. Med.* 329, 2002-2012 (1993).
5. Nathan, C. and Xie, Q.W. Regulation of biosynthesis of nitric oxide. *J. Biol. Chem.* 269, 13725-13728 (1994).
6. Schmidt, H.H., Lohmann, S.M., and Walter, U. The nitric oxide and cGMP signal transduction system: regulation and mechanism of action. *Biochim. Biophys. Acta* 1178, 153-175 (1993).
7. Stuehr, D.J., Cho, H.J., Kwon, N.S., *et al.* Purification and characterization of the cytokine-induced macrophage nitric oxide synthase: an FAD- and FMN-containing flavoprotein. *Proc. Natl. Acad. Sci. USA* 88, 7773-7777 (1991).
8. Schmidt, H.H., Pollock, J.S., Nakane, M., *et al.* Purification of a soluble isoform of guanylyl cyclase-activating-factor synthase. *Proc. Natl. Acad. Sci. USA* 88, 365-369 (1991).
9. Busconi, L. and Michel, T. Endothelial nitric oxide synthase. N-terminal myristoylation determines subcellular localization. *J. Biol. Chem.* 268, 8410-8413 (1993).
10. Pollock, J.S., Forstermann, U., Mitchell, J.A., *et al.* Purification and characterization of particulate endothelium-derived relaxing factor synthase from cultured and native bovine aortic endothelial cells. *Proc. Natl. Acad. Sci. USA* 88, 10480-10484 (1991).
11. Kobzik, L., Reid, M.B., Bredt, D.S., *et al.* Nitric oxide in skeletal muscle. *Nature* 372, 546-548 (1994).
12. Knowles, R.G., Palacios, M., Palmer, R.M., *et al.* Formation of nitric oxide from L-arginine in the central nervous system: a transduction mechanism for stimulation of the soluble guanylate cyclase. *Proc. Natl. Acad. Sci. USA* 86, 5159-5162 (1989).
13. Bredt, D.S. and Snyder, S.H. Isolation of nitric oxide synthetase, a calmodulin-requiring enzyme. *Proc. Natl. Acad. Sci. USA* 87, 682-685 (1990).
14. Salter, M., Knowles, R.G., and Moncada, S. Widespread tissue distribution, species distribution and changes in activity of Ca^{2+} -dependent and Ca^{2+} -independent nitric oxide synthases. *FEBS Lett.* 291, 145-149 (1991).

RELATED PRODUCTS

iNOS (murine macrophage) - Cat. No. 60862 • iNOS (murine recombinant) - Cat. No. 60864 • nNOS (rat recombinant) - Cat. No. 60870 • nNOS (rat recombinant - purified) - Cat. No. 60875 • eNOS (bovine recombinant) - Cat. No. 60880 • Nitrate/Nitrite Colorimetric Assay Kit (LDH Method) - Cat. No. 760871

QUICK REFERENCE PROTOCOL

Preparation of Extracts from Tissues and Cultured Cells

Extraction of Proteins from Tissues

- Add 20 volumes of ice-cold 1X homogenization buffer to a tissue sample.
- Homogenize the tissue and spin the homogenate in a microcentrifuge for 5 minutes at 4°C.
- Transfer the supernatant to a fresh tube and keep the tube on ice until use.

Extraction of Proteins from Tissue Culture Cells

- Remove the culture media from the tissue culture cells, wash the cells with PBS and harvest the cells in PBS containing 1 mM EDTA.
- Spin the cells at full speed for 2 minutes, remove the supernatant by vacuum aspiration and resuspend the cells in homogenization buffer.
- Sonicate or homogenize the cells.
- Spin the cells at full speed for 5 minutes, remove the supernatant, and adjust the protein sample to a concentration of 5-10 µg/ml.

Measurement of Nitric Oxide Synthase Activity in Enzyme Extracts

- Prepare the reaction mixture in a microcentrifuge tube and store on ice.
- Add the tissue extract to 40 µl of the reaction mixture.
- Incubate the reaction at 22-37°C for 10-60 minutes.
- Add 400 µl of stop buffer and 100 µl of equilibrated resin to the reaction sample.
- Transfer the reaction sample to a spin cup and spin at full speed for 30 seconds.
- Transfer the eluate to a scintillation vial, add scintillation fluid, and quantitate the radioactivity in a liquid scintillation counter.
- To determine the ratio of unreacted arginine to citrulline, add 400 µl of elution buffer to the spin cup, place the spin cup in a microcentrifuge tube, and spin at full speed for 30 seconds.
- Transfer the eluate to a scintillation vial, add scintillation fluid, and quantitate the radioactivity in a liquid scintillation counter.

NOTES

Poly-Prep® is a registered trademark of Bio-Rad Laboratories, Inc.

This document is copyrighted. All rights are reserved. This document may not, in whole or part, be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form without prior consent, in writing, from Cayman Chemical Company.

©03/15/2002 Cayman Chemical Company, Ann Arbor, MI, All rights reserved. Printed in U.S.A.